



## GreenLine®/D™ Fiber Optic Laryngoscopes Disposable Blades & LED Handles

### GreenLine®/D™ Green Heel Blades

**5 YEAR STERILITY**

- Reduce the risk of cross-contamination with single patient use blades
- Constructed of non-painted 303/304 surgical stainless steel
- Acrylic stem facilitates maximum LED illumination for superior, clear visualization
- Soft matte finish prevents reflection and glare during tissue visualization
- Durable nylon blade heel with two ball bearing connections
- Extended safety heel design prevents blade-to-handle contact
- Fully disposable; product recycling available



GreenLine®/D™ Blades



Not Made with Natural Rubber Latex

**STERILE EO**

### GreenLine®/D™ LED Handles

**5 YEAR STERILITY + SHELF LIFE**

- Reduce the risk of cross-contamination with single patient use handles
- Textured surface ensures secure grip
- Mercury-free, pre-installed batteries
- Contained unit is fully disposable; product recycling available



5-5338-11

5-5338-10

5-5338-09

5-5338-12



5-0236-69



5-0236-70



Not Made with Natural Rubber Latex

**STERILE EO**

Item	Size	Description	Sterile
5-5338-09	Medium	Textured Grip	Sterile
5-5338-10	Penlite	Textured Grip	Sterile
5-5338-11	Stubby	Textured Grip	Sterile
5-5338-12	Standard	Textured Grip	Sterile
5-0236-69	Large	Ribbed Grip	Non-Sterile
5-0236-70	Penlite	Smooth Grip	Non-Sterile



**SunOne™**  
SYSTEM

Single Patient Use,  
Fully Disposable  
Laryngoscopes




## GreenLine®/D™ Macintosh Blades

Item	Size	Length	Height	Description
5-5332-01	1	92mm	22mm	Infant
5-5332-02	2	100mm	21mm	Child
5-5332-03	3	130mm	22mm	Medium Adult
5-5332-35	3.5	144mm	22mm	Ext. Med. Adult
5-5332-04	4	155mm	25mm	Large Adult

## GreenLine®/D™ Miller Blades

Item	Size	Length	Height	Description
5-5333-20	00	67mm	10mm	Premature
5-5333-00	0	77mm	10mm	Neonate
5-5333-01	1	102mm	10mm	Infant
5-5333-02	2	153mm	12mm	Child
5-5333-03	3	195mm	14mm	Medium Adult
5-5333-04	4	205mm	18mm	Large Adult



## GreenLine®/D™ Robertshaw Blades

Item	Size	Length	Height	Description
5-5337-00	0	94mm	11mm	Neonate
5-5337-01	1	112mm	13mm	Infant

## GreenLine®/D™ Flex Tip Blades

Item	Size	Length	Height	Description
5-5382-03	3	195mm	18mm	Medium Adult
5-5382-04	4	205mm	18mm	Large Adult

**SunMed**

 ASK ABOUT PRODUCT RECYCLING

[www.SunOneSystem.com](http://www.SunOneSystem.com)