

# TUM-E-VAC® GASTRIC LAVAGE

LATEX FREE | SINGLE USE | DISPOSABLE | CE MARK

Tum-E-Vac® Gastric Lavage Kits are fully assembled, disposable kits that promote quick, efficient lavage of stomach contents. This completely closed system is designed to protect healthcare professionals from potential exposure to stomach contents during gastric lavage or decontamination procedures.



## TUM-E-VAC® GASTRIC LAVAGE KIT


- Pre-assembled kits are quick and easy to use, reducing the time for lavage procedures
- Closed system helps prevent exposure to stomach contents
- Available with activated charcoal suspension (with or without Sorbitol)
- Easy-to-fill top port
- Available with Universal Solution Spike (for use with pre-filled solution containers)
- No clean up or disinfecting procedures necessary
- All kits come standard with Bite Block

ITEM#	DESCRIPTION	PK
2054	Tum-E-Vac + 3 L Solution Bag + 32 Fr Lavage Tubing	10
2055	Tum-E-Vac + 3 L Solution Bag	10
2075	Tum-E-Vac + Universal Solution Spike	10
2060	Tum-E-Vac Bite Block	10

# LAVAGE TUBING

LATEX FREE | SINGLE USE | DISPOSABLE | CE MARK

Lavage Tubes feature an open bullet tip design with multiple eyelets to create free-flowing movement with low resistance.



## LAVAGE TUBING

- Soft clear tube with smooth tip
- Suction connector on proximal end
- Radiopaque stripe
- All tubing is 48" long with reference marks at 18", 20" and 22" from distal end
- Individually packaged

CAT#	DESCRIPTION	SIZE	PK
2124	48" Lavate Tubing	24 Fr	10
2128	48" Lavate Tubing	28 Fr	10
2132	48" Lavate Tubing	32 Fr	10
2136	48" Lavate Tubing	36 Fr	10