


MRI-SAFE BREATHING CIRCUIT

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

The MRI-safe transport circuit is lightweight and made with components safe for use in the MRI suite. Corrugated tubing is expandable to 9' for easier ventilation management during procedures.



MRI-SAFE BREATHING CIRCUIT

- Patient valve, PEEP valve and oxygen tubing constructed of MRI-safe materials
- 3 o'clock position of PEEP keeps valve out of patient's face
- For use with Smiths Medical Pneupac and other ventilators (*Consult ventilator manufacturer specifications prior to use*)
- Corrugating tubing expands generously to 9'
- Circuit available with or without PEEP valve

ITEM#	DESCRIPTION	PK
VC1000	MRI Transport Circuit with PEEP Valve	25
VC1000F	MRI Transport Circuit with PEEP Valve + Filter	25
VC1010	MRI Transport Circuit	25
VC1010F	MRI Transport Circuit + Filter	10
VP705	19 mm MRI Safe Adjustable PEEP Valve	20